

BLOOD DONOR STATION



HI-TECH MEDICAL EQUIPMENTS

Introduces a state of art of
"HI-CARE" Blood Donor Station.

The Hi-care Blood Donor Station is designed to assist and facilitate the work of medical staff. It can be adjusted electrically to form either a comfortable easy chair or a bed like coach. It requires less space than a conventional hospital bed. Hi-Care Blood Station is designed for the management of Blood Transfusion service by insuring safety, Comfort and Conveyance to the donor as well as Phlebotomist.

Blood Centres and transfusion services totally depend on voluntary donors, to ensure such a voluntary donor for physically and mentally well being, it is essentials for the blood centres and transfusion services to offer an atmosphere that is relaxed comforting and non-clinical.

Professional Blood Donor Station with Trendelenburg position during Vaso-Vegal attack by pressing single button.

COMPREHENSIVE FEATURES & BENEFITS :-

1. Three Powerful actuators facilitating easy movement of Back Rest, Seat tilt & Leg rest.
2. Micro- Processor based control box for fault proof operation .
3. Single Button operation for easy entry and exit
4. Tilt up Arm Rest provides flexibility for patient easy side entry
5. The Cushioned Foot rest can be adjusted to fit the length of any patient.
6. Durable high quality washable Upholstery provides comfort & Easy maintenance.
7. I.V. POLE attached with phlebotomist tray to keep accessories like cotton swabs, betadine and band-aid etc.
8. Blood collection tray is provided below the hand rest. It is having flexibility to shift from left to right and Vice-Versa.
9. Smooth & Lockable P.U. based Castors for silent move around without leaving floor marks.
10. Stable long Lasting metal frame design
11. Flexible arm " READING LAMP " to convenience medical Procedures.
12. Smooth & low noise operation will not disturb the donors rest.
13. An additional multipurpose table attached to the chair to save valuable space of centres.
14. Optional 360° Adjustable LCD TV with VGA Connector
15. Optional storage Cabinet.