

Tru-Cut® / Biopsy Needles

Tru-Cut® Biopsy Needles

Tru-Cut biopsy needles are designed for manual capture of high-quality tissue samples with minimal trauma to the patient.



Features include:

- 20mm specimen notch.
- Scalpel-sharp, surgical steel cutting edge retrieves cleanly cut, high-quality specimens.
- Removable stylet allows multiple sampling capability.
- Thin-wall needle yields larger core.
- Centimeter depth markings assist in needle placement.

TRU-CUT BIOPSY NEEDLES

Cat. No.	Size
2N2711X	18G x 3"
2N2712X	18G x 4 1/2"
2N2713X	18G x 6"
2N2700X	14G x 3 "
2N2702X	14G x 4 1/2"
2N2704X	14G x 6"

For single use only. 10 needles per case.

BREAST	KIDNEY	LIVER	LUNG	PROSTATE	THYROID
--------	---------------	--------------	------	----------	---------