




Dialysis Products



COVIDIEN

Quinton Catheter Accessories




Quinton Faller Tunneling Stylet

Code	Description	Ship Case
8888415679	Quinton Faller Tunneling Stylet	1



Quinton 2 Part Titanium Connector

Code	Description	Ship Case
8888415604	Quinton 2 Part Titanium Connector	1



Quinton Peritoneal Tunneling Stylet


- Made of pliable stainless steel
- Can be bent to create customized, curved tunnel
- Sterile

Code	Description	Ship Case
8817791001	Quinton 5Fr. x 23cm (9") stylet *	5



Quinton Titanium Catheter Extender

Code	Description	Ship Case
8888415612	Quinton Adult Titanium Catheter Extender	1



Quinton Peritoneal Catheter Straightening Stylet

- Stainless steel stylet
- Used to keep the catheter straight and maneuverable during insertion into the peritoneum
- Autoclavable
- Non-sterile

Code	Description	Ship Case
8888415646	Quinton Insertion/Straightening Stylet, 44.5cm	1
8888415661	Quinton Insertion/Straightening Stylet, 61cm	1