




Dialysis Products



COVIDIEN

Short-Term Vascular Access



Mahurkar Straight Extension Dual Lumen Catheters

- Made of flexible radiopaque polyurethane with a radiopaque tip
- Double-D catheter lumina
- Rotatable transparent suture wings
- Catheter extensions imprinted with name, length and priming volumes
- Two injection sealing caps
- Sterile

Code Side-holes	Code Side Slots	Catheter Size	Catheter Implant Length	Arterial Priming Volume	Venous Priming Volume	Ship Case
—	8813816005	11.5Fr.	13.5cm	1.0cc	1.1cc	5
—	8830414001	11.5Fr.	16cm	1.1cc	1.2cc	5
—	8813794005	11.5Fr.	19.5cm	1.2cc	1.3cc	5
—	8831662001	11.5Fr.	24cm	1.4cc	1.5cc	5
8817142005	—	10Fr.	12cm	0.8cc	0.9cc	5
8817145005	—	10Fr.	15cm	0.9cc	1.0cc	5



Mahurkar Straight Extension Dual Lumen Catheter Kits

Sterile Contents:

- (1) 8, 10 or 11.5Fr. Catheter
- (1) 18 ga. Introducer Needle
- (1) J/Straight Guidewire (0.032"/8Fr., 0.035"/10Fr. or 0.038"/11.5Fr.)
- (1) 8Fr. or 10Fr. Dilator
- (1) 12Fr. Dilator (11.5Fr. Kits Only)
- (2) Telfa Island Dressings
- (2) Injection Sealing Caps

Code Side-holes	Code Side Slots	Catheter Size	Catheter Implant Length	Arterial Priming Volume	Venous Priming Volume	Ship Case
—	8813817005	11.5Fr.	13.5cm	1.0cc	1.1cc	5
—	8830415001	11.5Fr.	16cm	1.1cc	1.2cc	5
—	8813793009	11.5Fr.	19.5cm	1.2cc	1.3cc	5
—	8831661001	11.5Fr.	24cm	1.4cc	1.5cc	5
8832539001	—	8Fr.	9cm	0.8cc	0.9cc	5
8832539002	—	8Fr.	12cm	0.9cc	0.9cc	5
8832539003	—	8Fr.	15cm	1.0cc	1.0cc	5
8817146001	—	10Fr.	15cm	0.9cc	1.0cc	5

Short-Term Vascular Access



Mahurkar Curved Extension Dual Lumen Catheters

- Patented curved extensions for jugular placement
- 180° curve bends away from patient's neck
- Made of flexible radiopaque polyurethane with a radiopaque tip
- Double-D catheter lumina
- Rotatable transparent suture wings
- Catheter extensions imprinted with name, length and priming volumes
- Two injection sealing caps
- Sterile

Code Side Slots	Catheter Size	Catheter Implant Length	Arterial Priming Volume	Venous Priming Volume	Ship Case
8813816009	11.5Fr.	13.5cm	1.1cc	1.2cc	5
8830414002	11.5Fr.	16cm	1.2cc	1.3cc	5
8813794009	11.5Fr.	19.5cm	1.3cc	1.4cc	5



Mahurkar Curved Extension Dual Lumen Catheter Kits

Sterile Contents:

- (1) 8, 10 or 11.5Fr. Catheter
- (1) 18 ga. Introducer Needle
- (1) J/Straight Guidewire (0.032", 0.035" or 0.038")
- (1) 8Fr. or 10Fr. Dilator
- (1) 12Fr. Dilator (11.5Fr. Kits Only)
- (2) Telfa Island Dressings
- (2) Injection Sealing Caps

Code Side-holes	Code Side Slots	Catheter Size	Catheter Implant Length	Arterial Priming Volume	Venous Priming Volume	Ship Case
—	8813817009	11.5Fr.	13.5cm	1.1cc	1.2cc	5
—	8830415003	11.5Fr.	16cm	1.2cc	1.3cc	5
—	8813793013	11.5Fr.	19.5cm	1.3cc	1.4cc	5
8832539006	—	8Fr.	12cm	1.0cc	1.1cc	5
8832539007	—	8Fr.	15cm	1.1cc	1.2cc	5
8817143005	—	10Fr.	12cm	0.9cc	1.0cc	5
8817146007	—	10Fr.	15cm	1.0cc	1.1cc	5
8817149007	—	10Fr.	19.5cm	1.1cc	1.2cc	5

Short-Term Vascular Access

Mahurkar 12Fr. Triple Lumen Catheters

- Soft Atraumatic Tip designed to reduce vessel trauma during insertion
- Laser-Cut side slots help to reduce positional occlusion
- Imprinted Priming Volumes on catheter extensions quickly identify heparin locking volumes
- Modified Double-D Lumen design ensures high flows
- 180° Curved Extensions help stabilize the catheter and provide patient comfort
- Rotatable Suture Wings enable the tip to be rotated after implantation and stabilizes the catheter



Mahurkar 12Fr. Triple Lumen Curved Extension Catheter Kits

Sterile Contents:

- (1) 12Fr. Catheter
- (1) 18 ga. Introducer Needle
- (1) J/Straight 0.038" Guidewire
- (1) 10Fr. Dilator
- (1) 12Fr. Dilator
- (2) Telfa Island Dressings
- (3) Injection Sealing Caps



Modified
Double-D
Lumen
Design

Code Side Slot	Catheter Size	Catheter Implant Length	Lumen Priming Volumes			Ship Case
			Blue or Distal	Red or Proximal	Clear or Medial	
8888345603	12Fr.	13cm	1.4cc	1.4cc	0.4cc	5
8888345611	12Fr.	16cm	1.5cc	1.5cc	0.4cc	5
8888345629	12Fr.	20cm	1.7cc	1.6cc	0.5cc	5
8888345637	12Fr.	24cm	1.8cc	1.8cc	0.5cc	5