

Jamshidi® Bone Marrow / Aspiration Needles

The original Jamshidi® Bone Marrow Biopsy/Aspiration Needles


A **lightweight handle** with built-in comfort knob enables the practitioner to apply sufficient pressure without discomfort. The needles are packaged with an **adapter to accommodate Luer lock syringes**.

The **tapered distal tip** facilitates the recovery of bone marrow samples with unaltered architecture.

The **single-use design** helps assure sterility for every procedure and helps prevent cross contamination. The needles are available in a **wide variety of sizes** to meet users' preference.

Ordering Information

Cat. No.	Description	Qty. Cs
DJ2013X	13G x 2"	10
DJ3513X	13G x 3 1/2"	10
DJ4011X	11G x 4"	10
DJ6011X	11G x 6"	10
DJ4008X	8G x 4"	10



Needle retrieves high-quality samples.



The razor-sharp cutting edge allows for easy coring of the bone.



Illinois Sternal/Iliac Bone Marrow/Aspiration Needles

The **sharp lancet point** penetrates easily into the bone to help ensure a safe and simple aspiration of marrow from the iliac crest or sternum.

An **adjustable depth guard** provides depth control for greater patient safety. The **twist-off cap** holds the stylet securely in place.

Our device is **compatible with Luer slip and Luer lock syringes** for a secure fit and ease in aspiration.

Ordering Information

Cat. No.	Description	Qty. Cs
DIN1515X	15G (adjustable length 3/8" to 17/8")	10
DIN1518X	15G (adjustable length 3/16" to 17/16")	10

