

COMPRESSOR NEBULIZER

For Better Health Care



Easy To Inhale

Hands Free Operation

Easy To Carry



ASPEN DIAGNOSTICS (P) LTD.
(Rapid Diagnostic Group of Companies)
ISO 9001 : 2000 CERTIFIED COMPANY

CE 0051

COMPRESSOR NEBULIZER

A-410

For Better Health Care

STANDARD ACCESSORIES PROVIDED:

1. Nebulizer (Medicine Cup)	1 Pc.
2. Mouth Piece	1 Pc.
3. Pediatric & Adult Mask	1 Pc. (each)
4. Nose Piece	1 Pc.
5. Air Connection Hose (Tube)	1 Pc.
6. Filters	2 Pcs.



SPECIFICATIONS

Model	Compressor Nebulizer A-410	Maximum Pressure:	2.5 (bar)
Dimensions:	Approx. 170 (W) x 150 (H) x 95 (D) mm	Maximum Air Flow:	10 L/min
Weight:	Approx. 1.7 kg	Medication Capacity:	8 ml
Sound Level:	55 dBa (with Nebulizer)	Nabulizing Rate:	0.35 ml / min.
Power Supply	230V ~ 50 Hz	Operating Temperature:	10°C to 40°C
Particle Size:	< 5 µm	Storage Temperature:	-25°C to +70°C

Import & Mktd. by:

ASPEN DIAGNOSTICS (P) LTD.
(Rapid Diagnostic Group of Companies)
ISO 9001 : 2000 CERTIFIED COMPANY

Mfd. by: E.P.R., Italy

