



Curity™ Silicone Coated Latex Foley Catheters



COVIDIEN

- Silicone coating provides a barrier between the latex catheter shaft and the patient's urethra increasing biocompatibility
- Large drainage eyes and drainage lumen facilitates urine flow
- Thick funnel helps prevent collapse when irrigating
- Strong symmetric balloons enhances confidence in catheter reliability
- For short term use

Description Size	Series 3			Series 1			Series 6			Series 7			Series 03			Series 04**		
	Hard Valve Paper Pouch			Hard Valve Poly Pouch			Rubber Valve Paper Pouch			Rubber Valve Poly Pouch			Rubber Valve Poly Pouch			Rubber Valve Poly Pouch		
Paediatric	Code	Qty per Carton	Qty per Case	Code	Qty per Bag	Qty per Case	Code	Qty per Bag	Qty per Case	Code	Qty per Bag	Qty per Case	Code	Qty per Carton	Qty per Case	Code	Qty per Bag	Qty per Case
8 Fr x 3cc	11083	12	1200	11081	10	800	11086	10	1000	11087	10	800	2908-03	10	1000	2908-04	100	2000
10 Fr x 3cc	11103	12	1200	11101	10	800	11106	10	1000	11107	10	800	2910-03	10	1000	2910-04	100	2000
2 Way	Code	Qty per Carton	Qty per Case	Code	Qty per Bag	Qty per Case	Code	Qty per Bag	Qty per Case	Code	Qty per Bag	Qty per Case	Code	Qty per Carton	Qty per Case	Code	Qty per Bag	Qty per Case
12 Fr x 5cc	25123	12	1200	25121	10	800	25126	10	1000	25127	10	800	1125-03	10	1000	1125-04	100	2000
14 Fr x 5cc	25133	12	1200	25131	10	800	25136	10	1000	25137	10	800	1185-03	10	1000	1185-04	100	2000
16 Fr x 5cc	25163	12	1200	25161	10	800	25166	10	1000	25167	10	800	1213-03	10	1000	1213-04	100	1600
18 Fr x 5cc	25183	12	1200	25181	10	700	25186	10	1000	25187	10	700	1394-03	10	1000	1394-04	100	1400
20 Fr x 5cc	25193	12	1200	25191	10	600	25196	10	1000	25197	10	600	1440-03	10	1000	1440-04	100	1200
22 Fr x 5cc	25223	12	1200	25221	10	600	25226	10	1000	25227	10	600	1575-03	10	1000	1575-04	100	1200
24 Fr x 5cc	25243	12	1200	25241	10	500	25246	10	1000	25247	10	500	1620-03	10	1000	1620-04	100	1200
26 Fr x 5cc	25263	12	1200	25261	10	500	25266	10	1000	25267	10	500	1746-03	10	1000	1746-04	100	1200
28 Fr x 5cc	25283	12	1200	25281	10	500	25286	10	1000	25287	10	500	1772-03	10	1000	1772-04	100	1200
30 Fr x 5cc	25303	12	1200	25301	10	500	25306	10	1000	25307	10	500	1800-03	10	1000	1800-04	100	1200
12 Fr x 30cc	26123	12	1200	26121	10	800	26126	10	1000	26127	10	800	5077-03	10	1000	5077-04	100	2000
14 Fr x 30cc	26143	12	1200	26141	10	800	26146	10	1000	26147	10	800	5101-03	10	1000	5101-04	100	2000
16 Fr x 30cc	26163	12	1200	26161	10	800	26176	10	1000	26167	10	800	5230-03	10	1000	5230-04	100	1600
18 Fr x 30cc	26183	12	1200	26181	10	700	26186	10	1000	26187	10	700	5353-03	10	1000	5353-04	100	1400
20 Fr x 30cc	26203	12	1200	26201	10	600	26206	10	1000	26207	10	600	5440-03	10	1000	5440-04	100	1200
22 Fr x 30cc	26223	12	1200	26221	10	600	26226	10	1000	26227	10	600	5588-03	10	1000	5588-04	100	1200
24 Fr x 30cc	26243	12	1200	26241	10	500	26246	10	1000	26247	10	500	5611-03	10	1000	5611-04	100	1200
26 Fr x 30cc	26263	12	1200	26261	10	500	26266	10	1000	26267	10	500	5729-03	10	1000	5726-04	100	1200
28 Fr x 30cc	26273	12	1200	26271	10	500	26276	10	1000	26277	10	500	5830-03	10	1000	5830-04	100	1200
30 Fr x 30cc	26293	12	1200	26291	10	500	26296	10	1000	26297	10	500	5909-03	10	1000	5909-04	100	1200
3 Way	Code	Qty per Carton	Qty per Case	Code	Qty per Bag	Qty per Case	Code	Qty per Bag	Qty per Case	Code	Qty per Bag	Qty per Case	Code	Qty per Carton	Qty per Case	Code	Qty per Bag	Qty per Case
16 Fr x 5cc	11163	12	1200	11161	10	800	11166	10	1000	11167	10	800	2939-03	10	1000	2939-04	100	1600
18 Fr x 5cc	11183	12	1200	11181	10	700	11186	10	1000	11187	10	700	2940-03	10	1000	2940-04	100	1400
20 Fr x 5cc	11203	12	1200	11201	10	600	11206	10	1000	11207	10	600	2941-03	10	1000	2941-04	100	1200
22 Fr x 5cc	11223	12	1200	11221	10	600	11226	10	1000	11227	10	600	2942-03	10	1000	2942-04	100	1200
24 Fr x 5cc	11243	12	1200	11241	10	500	11246	10	1000	11247	10	500	2943-03	10	1000	2943-04	100	1200
26 Fr x 5cc	11263	12	1200	11261	10	500	11266	10	1000	11267	10	500						
16 Fr x 30cc	12163	12	1200	12161	10	800	12166	10	1000	12167	10	800	9716-03	10	1000	9716-04	100	1600
18 Fr x 30cc	12183	12	1200	12181	10	700	12186	10	1000	12187	10	700	9718-03	10	1000	9718-04	100	1400
20 Fr x 30cc	12203	12	1200	12201	10	600	12206	10	1000	12207	10	600	9720-03	10	1000	9720-04	100	1200
22 Fr x 30cc	12223	12	1200	12221	10	600	12226	10	1000	12227	10	600	9722-03	10	1000	9722-04	100	1200
24 Fr x 30cc	12243	12	1200	12241	10	500	12246	10	1000	12247	10	500	9724-03	10	1000	9724-04	100	1200
26 Fr x 30cc	12263	12	1200	12261	10	500	12266	10	1000	12267	10	500	9726-03	10	1000	9726-04	100	1200

Note: Only sold in Case quantities

** Series 04 is unsterile

Specifications subject to change without notice.



COVIDIEN and COVIDIEN with logo are trademarks of Covidien AG.
© 2007 Covidien AG or its affiliates. All rights reserved.
SG-UR-CUR 001 10/2007.

TYGO HEALTHCARE PTE LTD.
NO.26, ANG MO KIO
INDUSTRIAL PARK 2, #04-01
SINGAPORE 569507

65.6482.0100 (T)
65.6482.0300 (B)

WWW.COVIDIEN.COM